

How to install Lasergene 17.x Client:

Windows:

1. Download DNASTARLasergene17xWinInstaller_AE.zip
2. Open the zip archive
3. Copy the LasergeneClient17 folder to desktop
4. Open the folder
5. Run the application DNASTARLasergene17xxWinInstaller.exe
6. Follow the steps in the installation dialog to install the Lasergene client

Mac:

1. Download DNASTARLasergene17xxMacInstaller_AE.dmg to desktop
2. Open the DMG file
3. Open the LasergeneClient17 folder
4. Run the application DNASTAR Lasergene Installer
5. Follow the steps in the installation dialog to install the software